



Test kit, nitrogen, nitrate color disc, model NI-14

Številka izdelka: 1416100-EU

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO₃ and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Color disc/ Cadmium reduction
Merilno območje:	Dual Range
Model:	NI-14
Najmanjši korak rasti:	0,05
Parameter:	Nitrogen, Nitrate - As NO ₃ ⁻ -N
Platforma:	Color Disc
Smallest increment:	0,05 0,5
Standardno merilno območje:	0 - 1 mg/L NO ₃ ⁻ -N; 0 - 10 mg/L NO ₃ ⁻ -N
Število testov:	50
Test Kit model:	NI-14**
Vrsta primera:	B