



Be Right™



## Test Kit, chlorine (free and total), model CN-80

Številka izdelka: 2129000-EU

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Simple, low-cost color disc method is the most accurate visual test kit method.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Specifikacije

Ime metode:	Color Disc/DPD
Merilno območje:	0 - 0.68 mg/L Cl <sub>2</sub>
Model:	CN-80
Parameter:	Chlorine, low range - As free & total Cl <sub>2</sub>
Platforma:	Color Disc
Smallest increment:	0.02
	0.1
	0.5
Standardno merilno območje:	0 - 0.68 mg/L Cl <sub>2</sub> ; 0 - 3.4 mg/L Cl <sub>2</sub> ; 0 - 10 mg/L Cl <sub>2</sub>
Število testov:	100
Test Kit model:	CN-80 <sup>2</sup>
Vrsta primera:	D