



Ammonia Ionic Strength Adjustor (ISA), Solution, 500 mL

Številka izdelka:

2824349-EU

Nearno

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan



Ionic Strength Adjustor (ISA) for Ammonia determinations by Ion Selective Electrode (ISE) method. Eliminate most interferences. Solution, 500 mL bottle.

Specifikacije

Koncentracija:	ISA
Package settings:	Bottle
Parameter:	Ammonium
Prostornina / velikost paketa:	500 mL