



Be Right™



Nitrate test kit, model NI-11

Številka izdelka: 146803-EU

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



One of Hach's most popular test kits, it uses the simple, low-cost color disc method.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikacije

Ime metode:	Color disc/ Cadmium reduction
Merilno območje:	Single Range
Model:	NI-11
Parameter:	Nitrogen, Nitrate - As NO_3^- -N
Platforma:	Color Disc
Smallest increment:	1.0
Standardno merilno območje:	0 - 40 mg/L NO_3^- -N
Število testov:	100
Test Kit model:	NI-11**
Vrsta primera:	B